


Anti-Slc22A6 Rabbit Polyclonal Antibody For Wb And Elisa - Q4U2R8

Item Number: CM0019377



Introduction

Rabbit polyclonal antibody against human SLC22A6 (OAT1) organic anion transporter. Suitable for Western blot and ELISA, with cross-reactivity to human and rat. Ideal for renal drug transporter and pharmacokinetic research.

[Learn More](#)

Parameter	Value
Product Name	SLC22A6 Rabbit pAb
Remarks/Alias	OAT1; PAHT; HOAT1; ROAT1; SLC22A6
Species	Human
Gene ID (Human)	9356
Gene ID	9356
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-135 of human SLC22A6 (NP_004781.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q4U2R8
Protein Weight	62kDa
Shipping	Ice bag