

Anti-Slc22A7 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y694

Item Number: CM0012998

Introduction

High-quality rabbit polyclonal antibody against human SLC22A7 (OAT2/NLT), suitable for WB and ELISA. Cross-reacts with mouse and rat samples. Ideal for drug transport and renal/hepatic function research.

[Learn More](#)

Parameter	Value
Product Name	SLC22A7 Rabbit pAb
Remarks/Alias	NLT; OAT2; hOAT11; SLC22A7
Species	Human
Gene ID (Human)	10864
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-142 of human SLC22A7 (NP_696961.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y694
Protein Weight	60kDa
Shipping	Ice bag