

Anti-Slc23A2 Polyclonal Antibody For Wb, Elisa - Q9Ugh3

Item Number: CM0031031

Introduction

Rabbit polyclonal antibody targeting human SLC23A2 (SVCT2), suitable for Western blot and ELISA. Detects human and mouse. Ideal for investigating ascorbate transport and sodium-dependent vitamin C uptake.

[Learn More](#)

Parameter	Value
Product Name	SLC23A2 Rabbit pAb
Remarks/Alias	NBTL1; SVCT2; YSPL2; SLC23A1
Species	Human
GeneID (Human)	9962
GeneID	9962
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 551-650 of human SLC23A2 (NP_005107.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UGH3
Protein Weight	70kDa
Shipping	Ice bag