

Anti-Slc24A4 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nff2

Item Number: CM0023738

Introduction

Rabbit polyclonal antibody targeting SLC24A4 (NCKX4), validated for WB and ELISA; reacts with human, mouse, rat. Ideal for studies on calcium transport, olfaction, and amelogenesis.

[Learn More](#)

Parameter	Value
Product Name	SLC24A4 Rabbit pAb
Remarks/Alias	A12A5; NCKX4; SHEP6; SLC24A2; SLC24A4
Species	Human
GeneID	123041
Immunogen	Synthetic peptide corresponding to amino acids 200-300 of human SLC24A4 (NP_705932.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NFF2
Protein Weight	69kDa
Shipping	Ice bag