

Anti-Slc24A5 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q71Rs6

Item Number: CM0014383

Introduction

Rabbit polyclonal antibody targeting human SLC24A5, the sodium/potassium/calcium exchanger 5, for IF/ICC and ELISA. Involved in melanosome ion transport and pigmentation, useful for albinism research.

[Learn More](#)

Parameter	Value
Product Name	SLC24A5 Rabbit pAb
Remarks/Alias	JSX; OCA6; NCKX5; SHEP4; SLC24A5
Species	Human
GeneID (Human)	283652
GeneID	283652
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-300 of human SLC24A5 (NP_995322.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q71RS6
Protein Weight	55kDa
Shipping	Ice bag