

Anti-Slc25A1 Polyclonal Antibody For Wb, Elisa - P53007

Item Number: CM0022308

Introduction

Rabbit polyclonal antibody targeting human SLC25A1 (mitochondrial citrate transporter). Validated for WB and ELISA; cross-reacts with mouse and rat. Suitable for metabolic and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	SLC25A1 Rabbit pAb
Remarks/Alias	SLC25A1; CTP; D2L2AD; SEA; SLC20A3; solute carrier family 25 member 1
Species	Human
GeneID (Human)	6576
GeneID	6576
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 14-124 of human SLC25A1 (NP_005975.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53007
Protein Weight	34kDa
Shipping	Ice bag