

Anti-Slc25A1 Polyclonal Antibody For Wb, If/Icc, Elisa - P53007

Item Number: CM0018257

Introduction

Polyclonal antibody against human SLC25A1, validated for WB, IF/ICC, ELISA. Detects endogenous levels of SLC25A1 in Human, Mouse, Rat. Ideal for mitochondrial citrate transport and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A1 Rabbit pAb
Remarks/Alias	CIC; CTP; SEA; CMS23; D2L2AD; SLC20A3; SLC25A1
Species	Human
GeneID (Human)	6576
GeneID	6576
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-311 of human SLC25A1 (NP_005975.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53007
Protein Weight	34kDa
Shipping	Ice bag