

Anti-Slc25A11 Polyclonal Antibody For Wb, Elisa - Q02978

Item Number: CM0021204

Introduction

Research-grade rabbit polyclonal antibody targeting SLC25A11 (mitochondrial 2-oxoglutarate/malate carrier). Validated for WB and ELISA, cross-reacts with human and mouse. Ideal for mitochondrial metabolism and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A11 Rabbit pAb
Remarks/Alias	OGC, PGL6, SLC20A4, SLC25A11
Species	Human
Gene ID (Human)	8402
Gene ID	8402
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC25A11 (NP_003553.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q02978
Protein Weight	34kDa
Shipping	Ice bag