

Anti-Slc25A11 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q02978

Item Number: CM0028836

Introduction

Rabbit polyclonal antibody against SLC25A11 (OGC). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Detects mitochondrial 2-oxoglutarate/malate carrier protein.

[Learn More](#)

Parameter	Value
Product Name	SLC25A11 Rabbit pAb
Remarks/Alias	OGC; PGL6; SLC20A4; SLC25A11
Species	Human
GeneID (Human)	8402
GeneID	8402
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC25A11 (NP_003553.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02978
Protein Weight	34kDa
Shipping	Ice bag