


Anti-Slc25A12 Polyclonal Antibody For Wb, Elisa - 075746

Item Number: CM0016736



Introduction

Rabbit polyclonal antibody against human SLC25A12, validated for WB and ELISA. Reacts with human, mouse, and rat. Suitable for studying mitochondrial aspartate/glutamate transport and metabolic shuttles.

[Learn More](#)

Parameter	Value
Product Name	SLC25A12 Rabbit pAb
Remarks/Alias	AGC1; DEE39; ARALAR; EIEE39; SLC25A12
Species	Human
Gene ID (Human)	8604
Gene ID	8604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human SLC25A12 (NP_003696.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075746
Protein Weight	75kDa
Shipping	Ice bag