


Anti-Slc25A12 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Ip, Elisa - O75746

Item Number: CM0002751



Introduction

Recombinant rabbit monoclonal antibody targeting human SLC25A12, a mitochondrial aspartate/glutamate antiporter. Validated for WB, IF-P, IHC-P, IP, and ELISA; cross-reacts with mouse and rat. Ideal for mitochondrial transport and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	SLC25A12 Rabbit mAb
Remarks/Alias	AGC1; DEE39; ARALAR; EIEE39; SLC25A12
Species	Human
Gene ID (Human)	8604
Gene ID	8604
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75746
Protein Weight	75 kDa
Shipping	Ice bag