

Anti-Slc25A13 Polyclonal Antibody For Wb, Elisa - Q9Ujs0

Item Number: CM0016336



Introduction

Anti-SLC25A13 polyclonal antibody validated for WB and ELISA in human, mouse, and rat samples. Suitable for studying mitochondrial aspartate/glutamate antiport and citrin function.

[Learn More](#)

Parameter	Value
Product Name	SLC25A13 Rabbit pAb
Remarks/Alias	CTLN2; NICCD; CITRIN; ARALAR2; SLC25A13
Species	Human
Gene ID (Human)	10165
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SLC25A13 (NP_001153682.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJS0
Protein Weight	74kDa
Shipping	Ice bag