

Anti-Slc25A13/Citrin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Ujs0

Item Number: CM0004175



Introduction

High-quality rabbit monoclonal antibody against human SLC25A13/citrin (SWISS: Q9UJS0). Validated for WB, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for mitochondrial transport and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A13/Citrin Rabbit mAb
Remarks/Alias	SLC25A13; ARALAR2; CITRIN; CTLN2; solute carrier family 25 member 13; SLC25A13/Citrin
Species	Human
GeneID	10165
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human SLC25A13/Citrin (NP_055066.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJS0
Protein Weight	74kDa
Shipping	Ice bag