

Anti-Slc25A14 Rabbit Polyclonal Antibody For Wb, Elisa - O95258

Item Number: CM0017397

Introduction

SLC25A14 (UCP5/BMCP1) rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Suitable for mitochondrial carrier and metabolic research.

[Learn More](#)

Parameter	Details
Product Name	SLC25A14 Rabbit pAb
Remarks/Alias	UCP5; BMCP1; SLC25A14
Species	Human
GeneID	9016
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-200 of human SLC25A14 (NP_001269124.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95258
Protein Weight	36kDa
Shipping	Ice bag