

Anti-Slc25A16 Rabbit Polyclonal Antibody For Wb And Elisa - P16260

Item Number: CM0021763

Introduction

Rabbit polyclonal antibody targeting SLC25A16, a mitochondrial coenzyme A transporter linked to Graves disease. Suitable for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for mitochondrial and autoimmune research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A16 Rabbit pAb
Remarks/Alias	GDA; GDC; ML7; hML7; HGT.1; D10S105E; SLC25A16
Species	Human
GeneID	8034
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human SLC25A16 (NP_689920.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16260
Protein Weight	36kDa
Shipping	Ice bag