

Anti-Slc25A17 Polyclonal Antibody For Wb, If/Icc, Elisa - O43808

Item Number: CM0014372

Introduction

Rabbit polyclonal antibody to human SLC25A17 (PMP34), validated for WB, IF/ICC, and ELISA applications. Detects SLC25A17 in human and mouse, a peroxisomal cofactor transporter. Suitable for peroxisome biology research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A17 Rabbit pAb
Remarks/Alias	PMP34; SLC25A17
Species	Human
GeneID	10478
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human SLC25A17 (NP_006349.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O43808
Protein Weight	35kDa
Shipping	Ice bag