

Anti-Slc25A20 Polyclonal Antibody For Wb, Elisa - 043772

Item Number: CM0013411

Introduction

Rabbit polyclonal antibody targeting human SLC25A20, the mitochondrial carnitine/acylcarnitine carrier (CACT), validated for WB and ELISA, with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	SLC25A20 Rabbit pAb
Remarks/Alias	CAC; CACT; SLC25A20
Species	Human
GeneID	788
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC25A20 (NP_000378.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43772
Protein Weight	33kDa
Shipping	Ice bag