

Anti-Slc25A20 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43772

Item Number: CM0028145

Introduction

Rabbit polyclonal antibody against SLC25A20 (CAC/CACT), validated for WB, IF/ICC, ELISA in human and rat. Targets the mitochondrial carnitine/acylcarnitine carrier, essential for fatty acid oxidation.

[Learn More](#)

Parameter	Value
Product Name	SLC25A20 Rabbit pAb
Remarks/Alias	CAC; CACT; SLC25A20
Species	Human
Gene ID (Human)	788
Gene ID	788
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-301 of human SLC25A20 (NP_000378.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O43772
Protein Weight	33kDa
Shipping	Ice bag