

Anti-Slc25A21 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bqt8

Item Number: CM0022132

Introduction

SLC25A21 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Detects human and rat mitochondrial 2-oxodicarboxylate carrier (ODC), 33kDa. Ideal for metabolic transport and lysine catabolism research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A21 Rabbit pAb
Remarks/Alias	ODC; ODC1; MTDPS18; SLC25A21
Species	Human
GeneID (Human)	89874
GeneID	89874
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 148-247 of human SLC25A21 (NP_085134.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BQT8
Protein Weight	33kDa
Shipping	Ice bag