

Anti-Slc25A24 Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Nuk1

Item Number: CM0016563



Introduction

SLC25A24 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Targets mitochondrial ATP-Mg/Pi carrier SCaMC-1 (53 kDa). Ideal for studies on mitochondrial nucleotide transport and calcium homeostasis.

[Learn More](#)

Parameter	Value
Product Name	SLC25A24 Rabbit pAb
Remarks/Alias	APC1; SCAMC1; SCAMC-1; SLC25A24
Species	Human
GeneID (Human)	29957
GeneID	29957
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human SLC25A24 (NP_037518.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6NUK1
Protein Weight	53kDa
Shipping	Ice bag