

## Anti-Slc25A26 Rabbit Polyclonal Antibody For Wb And Elisa - Q70Hw3

Item Number: CM0020841

### Introduction

Rabbit polyclonal antibody against human SLC25A26, suitable for WB and ELISA, cross-reactive with mouse. Targets mitochondrial S-adenosylmethionine carrier protein for methylation and transport studies. Molecular weight: 29 kDa.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | SLC25A26 Rabbit pAb                                                                                                   |
| Remarks/Alias    | SAMC; COXPD28; SLC25A26                                                                                               |
| Species          | Human                                                                                                                 |
| GeneID           | 115286                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 105-185 of human SLC25A26 (NP_775742.4) |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Mouse                                                                                                                 |
| SWISS            | Q70HW3                                                                                                                |
| Protein Weight   | 29kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |