

Anti-Slc25A27 Polyclonal Antibody For Wb, If/Icc, Elisa - O95847

Item Number: CM0013329

Introduction

Rabbit polyclonal antibody targeting human SLC25A27 (UCP4). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Useful for mitochondrial uncoupling and metabolic stress studies.

[Learn More](#)

Parameter	Value
Product Name	SLC25A27 Rabbit pAb
Remarks/Alias	UCP4; SLC25A27
Species	Human
GeneID (Human)	9481
GeneID	9481
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 201-300 of human SLC25A27 (NP_004268.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95847
Protein Weight	36kDa
Shipping	Ice bag