

Anti-Slc25A38 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Dw6

Item Number: CM0013515

Introduction

Polyclonal antibody targeting human SLC25A38, a mitochondrial glycine transporter involved in heme biosynthesis. Validated for WB, IF/ICC, ELISA with human and mouse samples. Ideal for erythropoiesis and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A38 Rabbit pAb
Remarks/Alias	SIDBA2; SLC25A38
Species	Human
GeneID (Human)	54977
GeneID	54977
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human SLC25A38 (NP_060345.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96DW6
Protein Weight	34kDa
Shipping	Ice bag