

Anti-Slc25A43 Polyclonal Antibody For Wb, Elisa - Q8Wut9

Item Number: CM0018624

Introduction

Anti-SLC25A43 rabbit polyclonal antibody for WB and ELISA applications. Recognizes human solute carrier family 25 member 43, a mitochondrial inner membrane transporter. Ideal for mitochondrial transport studies.

[Learn More](#)

Parameter	Specification
Product Name	slc25a43 Rabbit pAb
Remarks/Alias	SLC25A43; slc25a43
Species	Human
GeneID (Human)	203427
GeneID	203427
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human slc25a43 (NP_660348.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WUT9
Protein Weight	38kDa
Shipping	Ice bag