

Anti-Slc25A4/Ant1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P12235

Item Number: CM0003159

Introduction

Rabbit monoclonal antibody against human SLC25A4/ANT1 (ADP/ATP translocase 1), KO validated for WB, IHC-P, IF/ICC, IP, and ELISA. Detects human, mouse, and rat. Ideal for mitochondrial transport and bioenergetics research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SLC25A4/ANT1 Rabbit mAb
Remarks/Alias	T1; ANT; AAC1; ANT1; PEO2; PEO3; ANT 1; PEOA2; MTDPS12; MTDPS12A
Species	Human
GeneID (Human)	291
GeneID	291
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human SLC25A4/ANT1 (NP_001142.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12235
Protein Weight	33kDa
Shipping	Ice bag