

Anti-Slc25A4/Ant1 Polyclonal Antibody For Wb, If/Icc, Elisa - P12235

Item Number: CM0013982

Introduction

Rabbit polyclonal antibody targeting human SLC25A4/ANT1 (ADP/ATP translocase 1). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying ADP/ATP translocation, mitochondrial uncoupling, and related diseases.

[Learn More](#)

Parameter	Value
Product Name	SLC25A4/ANT1 Rabbit pAb
Remarks/Alias	T1; ANT; AAC1; ANT1; PEO2; PEO3; ANT 1; PEOA2; MTDPS12; MTDPS12A; SLC25A4/ANT1
Species	Human
GeneID (Human)	291
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12235
Protein Weight	33kDa
Shipping	Ice bag