

## Anti-Slc26A11 Polyclonal Antibody For Wb, Elisa - Q86Wa9

Item Number: CM0022314

### Introduction

Rabbit polyclonal antibody targeting human SLC26A11 (sodium-independent sulfate anion transporter) validated for WB and ELISA, cross-reacts with mouse. Useful for studies in kidney, brain, and transport biology.

[Learn More](#)

| Parameter        | Value                                                                                                                            |
|------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | SLC26A11 Rabbit pAb                                                                                                              |
| Remarks/Alias    | SLC26A11                                                                                                                         |
| Species          | Human                                                                                                                            |
| Gene ID          | 284129 (human)                                                                                                                   |
| Immunogen        | Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 500-600 of human SLC26A11 (NP_001159819.1) |
| Source           | Rabbit                                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                                            |
| Application      | WB, ELISA                                                                                                                        |
| Cross Reactivity | Mouse                                                                                                                            |
| SWISS            | Q86WA9                                                                                                                           |
| Protein Weight   | 65kDa                                                                                                                            |
| Shipping         | Ice bag                                                                                                                          |