

Anti-Slc26A2 Rabbit Polyclonal Antibody For Wb And Elisa - P50443

Item Number: CM0018261



Introduction

Rabbit polyclonal antibody targeting human SLC26A2, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying sulfate transport and cartilage development disorders.

[Learn More](#)

Parameter	Value
Product Name	SLC26A2 Rabbit pAb
Remarks/Alias	DTD; EDM4; DTDST; MST153; D5S1708; MSTP157; SLC26A2
Species	Human
GeneID (Human)	1836
GeneID	1836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human SLC26A2 (NP_000103.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50443
Protein Weight	82kDa
Shipping	Ice bag