


Anti-Slc26A3 Polyclonal Antibody For Wb And Elisa - P40879

Item Number: CM0028685



Introduction

Rabbit polyclonal antibody targeting SLC26A3 (DRA/CLD), a chloride-bicarbonate exchanger. Validated for Western blot and ELISA; useful for colon cancer and ion transport research.

[Learn More](#)

Parameter	Value
Product Name	SLC26A3 Rabbit pAb
Remarks/Alias	CLD, DRA, SLC26A3
Species	Human
Gene ID	1811
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 505-764 of human SLC26A3 (NP_000102.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P40879
Protein Weight	85kDa
Shipping	Ice bag