


Anti-Slc26A6 Polyclonal Antibody For Wb And Elisa - Q9Bxs9

Item Number: CM0018726



Introduction

Rabbit polyclonal antibody targeting human SLC26A6 (pendrin-like protein 1) for Western blot (WB) and ELISA. Validated for human samples, it supports research on anion exchange transporters in acid-base homeostasis and ion transport, especially in kidney and pancreas.

[Learn More](#)

Parameter	Value
Product Name	SLC26A6 Rabbit pAb
Remarks/Alias	SLC26A6
Species	Human
Gene ID (Human)	65010
Gene ID	65010
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 509-758 of human SLC26A6 (NP_599025.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BXS9
Protein Weight	83kDa
Shipping	Ice bag