


Anti-Slc26A9 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Lbe3

Item Number: CM0019285



Introduction

High-quality rabbit polyclonal antibody against human SLC26A9, an ion transporter and chloride channel. Validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for respiratory and epithelial transport research.

[Learn More](#)

Parameter	Value
Product Name	SLC26A9 Rabbit pAb
Remarks/Alias	SLC26A9
Species	Human
GeneID (human)	115019
GeneID	115019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-660 of human SLC26A9 (NP_443166.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7LBE3
Protein Weight	87kDa
Shipping	Ice bag