

Anti-Slc27A1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Pcb7

Item Number: CM0017050

Introduction

SLC27A1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA with human, mouse, rat cross reactivity. Investigate long-chain fatty acid transport protein 1 in lipid metabolism.

[Learn More](#)

Parameter	Value
Product Name	SLC27A1 Rabbit pAb
Remarks/Alias	FATP; FATP1; ACSVL5; FATP-1; SLC27A1
Species	Human
GeneID (Human)	376497
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-646 of human SLC27A1 (NP_940982.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6PCB7
Protein Weight	71kDa
Shipping	Ice bag