

Anti-Slc27A2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O14975

Item Number: CM0014083

Introduction

Rabbit polyclonal antibody against SLC27A2 (FATP2), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for fatty acid transport and metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	SLC27A2 Rabbit pAb
Remarks/Alias	VLCS; FATP2; VLACS; ACSVL1; FACVL1; hFACVL1; HsT17226; SLC27A2
Species	Human
Gene ID (Human)	11001
Gene ID	11001
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-573 of human FATP2 (NP_003636.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14975
Protein Weight	70kDa
Shipping	Ice bag