

Anti-Slc27A6 Polyclonal Antibody For Ihc-P, Elisa - Q9Y2P4

Item Number: CM0021963

Introduction

Rabbit polyclonal antibody to SLC27A6 / FATP6 for IHC-P and ELISA applications. Detects human and mouse SLC27A6, a long-chain fatty acid transporter and acyl-CoA ligase, 70 kDa. Ideal for heart, placental, and adrenal research.

[Learn More](#)

Parameter	Value
Product Name	SLC27A6 Rabbit pAb
Remarks/Alias	SLC27A6; ACSVL2; FACVL2; FATP6; VLCS-H1; solute carrier family 27 member 6
Species	Human
GeneID	28965
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 506-600 of human SLC27A6 (NP_054750.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2P4
Protein Weight	70kDa
Shipping	Ice bag