

Anti-Slc27A6 Polyclonal Antibody For Wb And Elisa - Q9Y2P4

Item Number: CM0026434

Introduction

Rabbit polyclonal antibody to SLC27A6 (FATP6), validated for WB and ELISA. Reacts with human and mouse. Used in fatty acid transport and metabolism research. UniProt Q9Y2P4.

[Learn More](#)

Parameter	Value
Product Name	SLC27A6 Rabbit pAb
Remarks/Alias	FATP6; ACSVL2; FACVL2; VLCS-H1; SLC27A6
Species	Human
GeneID (Human)	28965
GeneID	28965
Immunogen	Recombinant protein; Recombinant fusion protein containing amino acids 370-619 of human SLC27A6 (NP_054750.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y2P4
Protein Weight	70kDa
Shipping	Ice bag