

Anti-Slc28A2 Rabbit Polyclonal Antibody For Wb & Elisa - O43868

Item Number: CM0022982

Introduction

Rabbit polyclonal antibody against human SLC28A2 (CNT2), suitable for Western blot and ELISA. Validated in human, mouse, and rat. Ideal for purine nucleoside transport research.

[Learn More](#)

Parameter	Value
Product Name	SLC28A2 Rabbit pAb
Remarks/Alias	CNT2, HCNT2, SPNT1, HsT17153, SLC28A2
Species	Human
Gene ID (Human)	9153
Gene ID	9153
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 350-450 of human SLC28A2 (NP_004203.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43868
Protein Weight	72kDa
Shipping	Ice bag