


Anti-Slc28A3 Polyclonal Antibody For Wb, Elisa - Q9Has3

Item Number: CM0018177



Introduction

SLC28A3 Rabbit pAb is a polyclonal antibody validated for Western blot and ELISA, detecting human and mouse concentrative nucleoside transporter 3, a sodium-dependent nucleoside transporter involved in drug transport.

[Learn More](#)

Parameter	Value
Product Name	SLC28A3 Rabbit pAb
Remarks/Alias	CNT3; SLC28A3
Species	Human
Gene ID (Human)	64078
Gene ID	64078
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 612-691 of human SLC28A3 (NP_071410.1).)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HAS3
Protein Weight	77kDa
Shipping	Ice bag