


## Anti-Slc29A1 Rabbit Polyclonal Antibody For Wb, Elisa - Q99808

Item Number: CM0016086



### Introduction

Rabbit polyclonal antibody against human SLC29A1/ENT1, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for nucleoside transport and equilibrative transporter research.

[Learn More](#)

| Parameter        | Value                                                                                                      |
|------------------|------------------------------------------------------------------------------------------------------------|
| Product Name     | SLC29A1 Rabbit pAb                                                                                         |
| Remarks/Alias    | AUG; ENT1; SLC29A1                                                                                         |
| Species          | Human                                                                                                      |
| GeneID (Human)   | 2030                                                                                                       |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 200-300 of human SLC29A1 (NP_001071645.1) |
| Source           | Rabbit                                                                                                     |
| Category         | Polyclonal Antibodies                                                                                      |
| Application      | WB, ELISA                                                                                                  |
| Cross Reactivity | Mouse, Rat                                                                                                 |
| SWISS            | Q99808                                                                                                     |
| Protein Weight   | 50kDa                                                                                                      |
| Shipping         | Ice bag                                                                                                    |