

Anti-Slc29A3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bzd2

Item Number: CM0016793

Introduction

Rabbit polyclonal antibody against human SLC29A3/ENT3. Validated for WB and ELISA. Targets the nucleoside transporter involved in lysosomal and mitochondrial purine/pyrimidine salvage. For nucleoside metabolism and macrophage immunology research.

[Learn More](#)

Parameter	Value
Product Name	SLC29A3 Rabbit pAb
Remarks/Alias	ENT3; HJCD; PHID; HCLAP; SLC29A3
Species	Human
GeneID (Human)	55315
GeneID	55315
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human SLC29A3 (NP_060814.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BZD2
Protein Weight	52kDa
Shipping	Ice bag