

Anti-Slc2A13 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Qe2

Item Number: CM0030728

Introduction

Rabbit polyclonal antibody against SLC2A13 (HMIT), a proton myo-inositol cotransporter, validated in WB and ELISA. Reacts with human, mouse, rat. For neuroscience and transporter research.

[Learn More](#)

Parameter	Value
Product Name	SLC2A13 Rabbit pAb
Remarks/Alias	HMIT; SLC2A13
Species	Human
GeneID (Human)	114134
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 564-648 of human SLC2A13 (NP_443117.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96QE2
Protein Weight	70kDa
Shipping	Ice bag