

Anti-Slc2A4Rg Polyclonal Antibody For Wb, Elisa - Q9Nr83

Item Number: CM0022053



Introduction

High-quality rabbit anti-SLC2A4RG polyclonal antibody for Western blotting and ELISA, with cross-reactivity to mouse. Ideal for studying GLUT4 and Huntington disease gene regulation mechanisms.

[Learn More](#)

Parameter	Value
Product Name	SLC2A4RG Rabbit pAb
Remarks/Alias	GEF; HDBP1; HDBP-1; Si-1-2; Si-1-2-19; SLC2A4RG
Species	Human
GeneID (Human)	56731
GeneID	56731
Immunogen	A synthetic peptide corresponding to amino acids 100-200 of human SLC2A4RG (NP_064446.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NR83
Protein Weight	41kDa
Shipping	Ice bag