

Anti-Slc2A5 (Glut5) Rabbit Polyclonal Antibody For Wb, Elisa - P22732

Item Number: CM0015643

Introduction

Rabbit polyclonal antibody against SLC2A5 (GLUT5), validated for Western blotting (WB) and ELISA, with cross-reactivity to human and mouse. Suitable for fructose transporter research and metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	SLC2A5 Rabbit pAb
Remarks/Alias	GLUT5; GLUT-5; SLC2A5
Species	Human
GeneID (human)	6518
GeneID	6518
Immunogen	Synthetic peptide, a synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC2A5 (NP_003030.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P22732
Protein Weight	55kDa
Shipping	Ice bag