


Anti-Slc2A9 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nrm0

Item Number: CM0017306



Introduction

Rabbit polyclonal antibody targeting SLC2A9 (GLUT9/URATv1), a high-capacity urate transporter. Suitable for WB and ELISA applications, with cross-reactivity to mouse and rat. Ideal for urate transport, renal function, and metabolic disorder research.

[Learn More](#)

Parameter	Value
Product Name	SLC2A9 Rabbit pAb
Remarks/Alias	GLUT9; GLUTX; UAQTL2; URATv1; SLC2A9
Species	Human
GeneID (human)	56606
GeneID	56606
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 429-529 of human SLC2A9 (NP_064425.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NRM0
Protein Weight	59kDa
Shipping	Ice bag