

Anti-Slc30A1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y6M5

Item Number: CM0025046



Introduction

Rabbit polyclonal antibody against SLC30A1/Zinc transporter 1, for WB, IHC-P, ELISA in human, mouse, rat. Investigates zinc/calcium transport, cuproptosis, and heavy metal influx protection. Optimized for research use.

[Learn More](#)

Parameter	Value
Product Name	SLC30A1 Rabbit pAb
Remarks/Alias	ZNT1; ZRC1; SLC30A1
Species	Human
Gene ID	7779
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 328-507 of human SLC30A1 (NP_067017.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6M5
Protein Weight	55kDa
Shipping	Ice bag