

Anti-Slc30A6 (Zinc Transporter 6) Rabbit Polyclonal Antibody For Wb And Elisa - Q6Nxt4

Item Number: CM0017952

Introduction

Rabbit polyclonal antibody targeting human SLC30A6 (ZNT6), validated for Western blot (WB) and ELISA. Detects human, mouse, and rat zinc transporter 6 with a predicted molecular weight of 51 kDa. Suitable for research on zinc homeostasis and the secretory pathway.

[Learn More](#)

Parameter	Value
Product Name	SLC30A6 Rabbit pAb
Remarks/Alias	ZNT6; MST103; MSTP103; SLC30A6
Species	Human
GeneID (Human)	55676
GeneID	55676
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 252-461 of human SLC30A6 (NP_060434.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6NXT4
Protein Weight	51kDa
Shipping	Ice bag