

Anti-Slc31A1 Polyclonal Antibody For Wb, Elisa - O15431

Item Number: CM0017400



Introduction

Rabbit polyclonal antibody targeting human SLC31A1 (CTR1, copper transporter 1). Validated for WB and ELISA with cross-reactivity to human and mouse. Ideal for copper metabolism research.

[Learn More](#)

Parameter	Value
Product Name	SLC31A1 Rabbit pAb
Remarks/Alias	CTR1; NSCT; COPT1; SLC31A1
Species	Human
Gene ID (Human)	1317
Gene ID	1317
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human SLC31A1 (NP_001850.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15431
Protein Weight	21 kDa
Shipping	Ice bag