

Anti-Slc31A2 Polyclonal Antibody For Wb, Elisa - O15432

Item Number: CM0020058

Introduction

Rabbit polyclonal antibody against mouse Slc31a2 (copper transporter CTR2). Validated for Western blot (WB) and ELISA. Ideal for copper transport research. Human SWISS: O15432.

[Learn More](#)

Parameter	Value
Product Name	Slc31a2 Rabbit pAb
Remarks/Alias	CTR2; Slc31a2
Species	Mouse
GeneID (human)	1318
GeneID	20530
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse Slc31a2 (NP_001277447.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O15432
Protein Weight	16kDa
Shipping	Ice bag