

Anti-Slc32A1 Rabbit Polyclonal Antibody For Wb - Q9H598

Item Number: CM0028051

Introduction

Rabbit polyclonal antibody against human SLC32A1 (VGAT/VIAAT), suitable for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying inhibitory neurotransmission and synaptic vesicle loading.

[Learn More](#)

Parameter	Value
Product Name	SLC32A1 Rabbit pAb
Remarks/Alias	VGAT, VIAAT, SLC32A1
Species	Human
GeneID (Human)	140679
GeneID	140679
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-525 of human SLC32A1 (Q9H598)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H598
Protein Weight	57kDa
Shipping	Ice bag