

Anti-Slc32A1/Vgat Polyclonal Antibody For Wb And Elisa - Q9H598

Item Number: CM0014578

Introduction

SLC32A1/VGAT Rabbit Polyclonal Antibody validated for WB and ELISA, recognizing human, mouse, and rat vesicular inhibitory amino acid transporter. Suitable for inhibitory neurotransmission research.

[Learn More](#)

Parameter	Value
Product Name	SLC32A1 Rabbit pAb
Remarks/Alias	VGAT; VIAAT; SLC32A1
Species	Human
GeneID (Human)	140679
GeneID	140679
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human SLC32A1 (NP_542119.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H598
Protein Weight	57kDa
Shipping	Ice bag