

Anti-Slc34A1 Rabbit Polyclonal Antibody For Wb And Elisa - Q06495

Item Number: CM0018377

Introduction

Rabbit polyclonal antibody against SLC34A1 (NaPi-2a), validated for WB and ELISA. Detects human, mouse, and rat renal sodium-phosphate cotransporter. Ideal for phosphate transport and kidney disease research.

[Learn More](#)

Parameter	Value
Product Name	SLC34A1 Rabbit pAb
Remarks/Alias	NPT2; FRTS2; SLC11; HCINF2; NAPI-3; NPTIIa; NPHLOP1; SLC17A2; SLC34A1
Species	Human
Gene ID (Human)	6569
Gene ID	6569
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 562-639 of human SLC34A1 (NP_003043.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06495
Protein Weight	69kDa
Shipping	Ice bag